



HERDEGEN "PROGRESS II" FULL CUFF CRUTCH

Crutches with adjustable elbow support especially designed to suit most user size and weight requirements. Easily adjustable with a push button, fully pivoting arm cuff. Adjustable elbow support height: 8.5 - 11.25 in (22-29 cm). Height adjustment between handle and floor: 28.25 - 37.75 in (72-96 cm). Lightweight 18.7 oz (530g) per crutch. Weight capacity: 330 lbs (150 kg) per crutch. Colour: white. For replacement crutch tips order 14A301-1.

Béquille canadienne « Progress II » de Herdegen

Une « béquille canadienne » avec un appui d'avant-bras réglable, conçue pour répondre à la plupart des exigences de taille et de poids d'un utilisateur. Muni d'un brassard pivotant facilement réglable en hauteur de 8,5 à 11,25 po (22 à 29 cm) grâce aux boutons-poussoirs. La béquille est ajustable en hauteur du sol à la poignée de 28,25 à 37,75 po (72 à 96 cm). Léger — chaque béquille pèse 18,7 oz (530 g). Capacité maximale de 330 lb (150 kg) par béquille. Couleur : blanc. Pour les embouts de rechange, commander 14A301-1.



14A301A



Limited One Year Warranty

Parsons ADL warrants to the original owner of this product that it is free from defects in material and workmanship for one year from date of purchase. This warranty is offered provided the product has not been abused, misused, altered or damaged after purchase. (This includes damage due to the use of harsh chemicals.)

Garantie limitée d'un an

Parsons ADL garantit au propriétaire d'origine du présent produit que celui-ci est exempt de défauts de matériaux et de fabrication pour une durée d'une année à compter de la date d'achat. Cette garantie est offerte à condition que le produit n'ait pas été utilisé abusivement, mal utilisé, modifié ou endommagé après l'achat. (Ceci comprend les dommages résultant de l'utilisation de produits chimiques puissants.)

Parsons ADL Inc. 1986 Sideroad 15, Tottenham, Ontario L0G 1W0 Canada

Tel: 905.936.3580 - 1.800.263-1281 / Fax: 905.936.3585 - 1.888.225.5535

Email: custserv@parsonsadl.com

Aids to Daily Living - Accessoires à la vie quotidienne