



Etac Turner Pro

Transfer platform - user turning with standing support

A sit to stand aid with functional design

Etac Turner Pro is a turn aid with a functional design that offers safe patient turning with standing support. This turn aid enables a good posture and minimises the risk of injury for the carer when transferring, moreover, it encourages natural standing and sitting movement for the user. The only adjustment necessary is the height of the leg support and the handle. Thereafter, the transfer can be performed safely whilst the caregiver and user can maintain eye contact and communication.

The orange handle provides a strong, high contrasting colour which can help users with impaired vision or dementia. The small footplate makes it the perfect choice of transfer aid when transferring in tight spaces. Etac Turner Pro is easy to transport due to its low weight of only 7.5 kg, furthermore it is compact and easy to grip.



Description

Etac Turner Pro

Dimensions

800 - 1200 mm

Item no.

16090106

Material

Cast aluminium, powder coated steel, chrome 3

Weight

Total weight 165 lbs
Maximum user weight (SWL) 440 lbs

Dimensions

Lower leg support height 150 mm
Lower leg support width 380 mm

Facts



Features

- The sit to stand aid is easy to transport due to a low weight of only 7.5 kg. It is compact, evenly balanced and easy to grip.
- The angled leg support is soft and comfortable and provides lateral support. Easy to adjust to whatever height is needed.
- The orange handle provides a strong, high contrasting colour which can help users with impaired vision or dementia.
- The cast aluminium base plate with integrated anti slip function is easy to turn. Together with the wheel, Turner Pro works well on different flooring materials including carpets.



“Everyone should be able to live a free and independent life
and pursue their dreams regardless of any physical circumstances”

For the latest news and continuously updated
product information – visit: www.etac.com

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 – 58 73 30
Fax +46 371 – 58 73 90
info@etac.se www.etac.com

 **etac**®
Creating Possibilities