



MULTIFOLD



- Separates into two pieces for easy carrying
- Rotation of pin either unlocks or locks the ramp sides together without extra tools
- Folds, carries like a suitcase
- Patented hinge eliminates pinch point
- Closure strap locks ramp panels together
- Lightweight, easy-to-handle & set up
- Durable welded construction
- Safety guide indicates if ramp is on a safe slope
- Full platform provides excellent stability and accommodates wheelchairs and scooters with various wheel configurations
- Anti-slip, high traction surface
- Safety DVD and steel security pins included
- Fits in the trunk of most cars for easy transportation



Item Number	Length	Width	Folded Size	Weight	Capacity*
WCR530	5'	30"	32¼" x 16" x 9"	30 lbs.	800 lbs.
WCR630	6'	30"	38¼" x 16" x 9"	35 lbs.	800 lbs.
WCR730	7'	30"	44¼" x 16" x 9"	39 lbs.	800 lbs.
WCR830	8'	30"	50¾" x 16" x 9½"	52 lbs.	800 lbs.
BAR536	5'	36"	32¼" x 17½" x 11½"	50 lbs.	800 lbs.
BAR636	6'	36"	38¼" x 17½" x 11½"	58 lbs.	800 lbs.
BAR736	7'	36"	44¼" x 17½" x 11½"	67 lbs.	800 lbs.
BAR836	8'	36"	50¾" x 17½" x 11½"	76 lbs.	800 lbs.


Not sure which ramp length you need? Refer to pages 2 - 3 to help you.



* Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles).



MULTIFOLD REACH



- Extended hook to clear rear bumper
- Red straps included for tailgate usage
- Separates into two pieces for easy carrying
- Rotation of pin either locks or allows ramp to be separated into two pieces
- Folds, carries like a suitcase
- Patented hinge eliminates pinch point
- Closure strap locks ramp panels together
- Lightweight, easy-to-handle & set up
- Durable welded construction
- Full platform provides excellent stability
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Safety DVD and steel security pins included
- Fits in the trunk of most cars for easy transportation



Safety Guide Bubble Level

Safety guide indicates if ramp is on a safe slope!

Item Number	Length	Width	Folded Size	Weight	Capacity*
UTW630	6'	30"	38¼" x 16" x 10¾"	36 lbs.	800 lbs.
UTW730	7'	30"	44¼" x 16" x 10¾"	41 lbs.	800 lbs.
UTW830	8'	30"	50¾" x 16" x 11"	52 lbs.	800 lbs.

Not sure which ramp length you need? Refer to pages 2 - 3 to help you.

* Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles).

