

Made in  
Germany

SPECIAL  
RECLINE  
BACK  
SUPPORT

## Aquatec *SRB*

Model no. 1573878

AQUATEC SRB, Special Reclining Backrest - White

The Aquatec *SRB* is equipped with adjustable laterals and padded chest strap to accommodate users with limited upper body control. The *SRB* also offers the ability to recline to 40 degrees and industry-leading 2 ½" minimum seat height, which makes it one of the most versatile bath lifts on the market.

## Features

- Compact construction with minimum seat height of 2 ½"
- Reinforced lifting mechanism
- Can be used with or without cover materials
- Ergonomic floating hand control
- Lithium-ion rechargeable battery pack
- Lift function will not allow the lift to descend if there is not enough battery power present to raise the user
- Can be easily folded and disassembled into two parts with one hand
- Self-releasing suction cups
- 40 degrees of recline

## Specifications

Seat Depth	19"
Seat Width	15"
Seat Width with Side Flaps Down	28"
Overall Length	24", 33" reclined
Back Height	26"
Back Width	16" - 27"
Base Width	11"
Base Length	23"
Minimum Seat Height	2 ½"
Maximum Lift Height (Optional Risers up to 19 ¾")	17"
Weight Capacity (lbs.)	300