



CHERRY RED



PARADISE BLUE



SOLAR FLARE ORANGE
Limited Edition



ENVY GREEN
Limited Edition



SATIN SILVER
Limited Edition



SUNBURST YELLOW
Limited Edition

Seat Color



Black Vinyl

**Dynamic LiNX Controller
Now Standard, Plus a New
Premium Seat Option!**

*Available in 6 Vibrant Colors
with a Larger Foot Plate
and New
Shroud Design*

**Disassembles
Quickly &
Easily!**

**LiteRider^{GOLDEN}
ENVY**
GP162


The all-new LiteRider Envy Power Wheelchair is both stylish and durable with a tubular design that is unique to the Golden LiteRider family! A Dynamic LiNX controller comes standard and the traditional style makes driving simple and reliable. The new premium seat offers contoured support with an increased back height of 17" for a longer seating tolerance and added comfort. The Envy, featuring 2-22 amp hour batteries for an operating range of up to 15.5 miles, is Golden's Retail Cash Power Wheelchair! It has a larger footplate, an impressive turning radius of 28.5", and disassembles quickly and easily for transport. When you are looking for a compact power wheel chair that is fun to drive and easy to transport – look no further than the Golden LiteRider Envy!



Model Number	GP162	
Medicare Code	Not Coded	
Weight Capacity	300 lbs.	
Drive Wheels	Rear	
Maximum Speed	3.5 mph	
Operating Range ¹	15.5 miles	
Ground Clearance	2.5"	
Turning Radius	28.5"	
Type Batteries	2-22 AH	
Colors	Cherry Red, Paradise Blue, Solar Flare Orange, Envy Green, Satin Silver, Sunburst Yellow	
Freewheel Mode	Yes	
Electronic Speed Control	Yes	
Electro-Mechanical Brakes	Yes	
Charger	Off Board	
Controller	40 Amp Dynamic Linx	
Dimensions:		
Length	36"	
Width	22.5"	
Height (ground to top of back of seat)	33" - 34"	
Height (ground to top of back of headrest)	N/A	
Ground to Top of Seat	20.5" - 21.5"	
Deck to Top of Seat	N/A	
Number of Seat Height Positions	2	
Size of Increments	1"	
Weight of Unit: (assembled)		
Front Half	123 lbs.	
Rear Half	37 lbs.	
Seat w/Arms	36 lbs.	
Seat w/Arms	23 lbs.	
Batteries	27 lbs/pack	
Tires:		
Drive Wheels	Flat Free	
Drive Wheels	9"	
Casters	6"	
Anti-Tip Wheels	3"	
Standard Seat:		
Back Height - no headrest	Stadium Style	
Back Height with headrest	16"	
Back Height with headrest	N/A	
Width x Depth	17" x 17"	
Color	Black Vinyl	
Optional Seats:		
Back Height - no headrest	Stadium Style	Premium
Back Height with headrest	16"	17"
Back Height with headrest	N/A	N/A
Width x Depth	20" x 18"	17" x 16"
Color	Black Vinyl	Black Vinyl
Seat Options:		
Pneumatic Rise Seat Option	N/A	
Electric Rise Seat Option	N/A	
Footrest:		
Height adjustable	Yes	
Angle adjustable	Yes	
Warranty:		
Frame	5 Years	
Drive Train	13 Months	
Electronics ²	13 Months	
Batteries	12 Months	

¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

² Electronics warranty excludes batteries.



Watch our video to learn about the benefits of the LiNX controller!

Scan here with smartphone camera.



Disassembles easily for transport!



Ergonomic LiNX controller



NEW! This model shown with Premium Seat Option!

Available Accessories:

- Cane Holder Tube
- Oxygen Tank Holder
- Cup Holder
- Quad Cane Holder
- Walker Holder
- Crutch Holder
- Safety Flag
- Power Chair Cover
- Pack N' Go
- Seat Belt (Optional)



Perfect for transport!

Some accessories require a mounting bracket or mounting clips.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

