

Please mark clearly with an X the items you wish to purchase.
Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744.
Contact our customer service team at (877) 252-6655.

Account #: _____	Ship to: _____
Date: _____	Name: _____
PO #: _____	Address: _____
Buyer: _____	City: _____
Attention: _____	Phone: () _____

Jogger Features:

- Super-light aluminium frame folds easily
- Comfortable padded seat reclines
- Adjustable 5-point padded chest harness and shoulder pads as well as a padded 3-point locking pelvic belt
- 12" pneumatic wheels with front swivel lock-out
- Removable machine-washable fabric



Jogger



Liners/Sitters

One-step, compact folding design
(Canopy folds with Jogger)



Front wheel rotates 360° and includes lock-out pin for rough terrain



Optional Accessories:



Optional rain cover



Soft-Touch liners and padded headrest



Soft-Touch Sitter



Removable canopy with window and height adjustable handle

Please mark clearly with an X the items you wish to purchase.
Email form to NuVision Rehab Group at sales@nuvisionrehab.com or fax form to 403-253-1744.
Contact our customer service team at (877) 252-6655.

Jogger Specifications	Stroller Only	Liner Size 1	Liner Size 2	Sitter Size 1	Sitter Size 2
Seat to Footwell	14"	16"	16"	17"	18"
Seat	14"W x 12"D				
Back Height	25"				
Maximum User Weight	110 lbs.				
Stroller Weight	27 lbs.				
Handle Height to Floor	35" - 43"				
Recline Range	100° - 140°				
Dimensions when Open	43"H x 26"W x 42"D				
Dimensions when Folded	16"H x 21"W x 33"D				

Item

Jogger

Part

62000000

Optional Accessories

Medical Necessities Bag (10 lb. weight capacity)

Part

62010000

Head Support

62020000

Rain Cover

62030000

Soft-Touch Liners Back Cushions

Size

Cherry 

Lilac 

Aqua 

Dark Gray 

LB1

72050101

72050102

72050104

72050109

LB2

72050201

72050202

72050204

72050209

Soft-Touch Liners Seat Cushions

Size

Cherry 

Lilac 

Aqua 

Dark Gray 

LS1

72040101

72040102

72040104

72040109

LS2

72040201

72040202

72040204

72040209

Soft-Touch Sitters

Size

Lilac 

Aqua 

Dark Gray 

Size 1

77010102

77010104

77010109

Size 2

77010202

77010204

77010209