



FEATURES

- Accommodate users weighing 20 - 102 lb and up to 60" tall
- Durable, machine washable, comfortable padded cover in three color options
- Crash tested to FMVSS #213
- Certified for use in aircraft
- Positioning options available
- Optional seat depth extender allows car seat to grow with the child

INCLUDES

- Height adjustable soft-padded headwings
- Memory foam back/seat padding reduces sores and fatigue
- Standard 5-point harness restraint for better protection than vehicle lap/shoulder belt
- Crotch strap pad
- Lateral positioning pads
- Dual-hook tether strap

UPHOLSTERY



Convertible Pink
(RR)



Roadster Red
(RR)



Speedway Gray
(SG)

ACCESSORIES



Extensor Thrust
Wedge



Trunk Support Pads



Buckle Guard



Retainer Clip Guard


ORDER FORM FAX TO 310-305-1718 OR SUBMIT ONLINE

Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____



BASE UNIT

CSI-2000P (Head pads)
\$1,195

CSI-2000W (Headwings)
\$1,195

Weight Capacity	20 -102 lb	20 -102 lb
Height Capacity	60"	60"
Shoulder Harness Heights	11.5", 14.5", 17.5", 20.5"	11.5", 14.5", 17.5", 20.5"
Back Height	29"	29"
Seat Depth	10"	10"
Seat Depth with Extender	12.5" or 15"	12.5" or 15"
Inside Seat Width at Front	15.5"	15.5"
Inside Seat Width at Rear	13"	13"
Overall Width	19"	19"
Overall Depth	11"	11"
Seat-to-Back Angle	88°	88°
Product Weight	11 lb	11 lb

UPHOLSTERY



Convertible Pink with Gray Mesh

Dune Buggy Beige with Gray Mesh

Go Kart Green with Gray Mesh

Roadster Red with Black Mesh

Speedway Gray with Gray Mesh



ACCESSORIES (fit both car seats)

Seat Depth Extender	CSI-A2025	Extend to 12.5"
	CSI-A2026	Extend to 15"
Abductor	CSI-A2015	For standard seat depth
	CSI-A2018	For 12.5" seat depth
	CSI-A2019	For 15" seat depth
Extensor Thrust Wedge	CSI-A2050	2" Wedge
	CSI-A2055	3" Wedge
Block Style Trunk Support Pads (2 each of 1.5" and 2.5" thick pads)	CSI-A2008	
Lateral Pads (set of 4)	CSI-A2010	
Bolt-in Hardware Kit with Strap(for vehicle manufactured before 2002)	CSI-A2031	
Buckle Guard	ZZ-A2075	
Retainer Clip Guard	CSZ-A2085	