

Platinum Collection Deluxe Padded Commode Bath Seats

Platinum Collection Deluxe Padded Commode Bath Seats

- Open front seat provides easier access for personal care
- Comfortable, cushioned seat and backrest are water-tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- Comes complete with commode pail and cover
- Available with or without support arms
- Comes assembled with tool-free arm and back attachments
- Maximum Weight Capacity 300 lb EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0245
- · Lifetime limited warranty

GRAHAM-FIELD

7946A without support arms 1 ea **7947A** with support arms 1 ea

<u>Specifications</u>	<u>7946A</u>	<u>7947A</u>
Overall Width	20 ½"	20 ½"
Overall Depth	21"	21"
Seat Width	15 ¾"	15 ¾"
Seat Depth	16"	16"
Minimum Seat Height	17 ½"	17 ½"
Maximum Seat Height	22 ½"	22 ½"
Height of Back Support Above Seat	16"	16"
Diameter Commode Opening	8"	8"
Width Between Arms	NA	17 ¾"
Maximum Weight Capacity 300 lb EV	ENLY DIS	STRIBUTED

