



SUITCASE[®]

Anytime. Anywhere.
Portable ramps for your every need.



EZ-ACCESS[®] SUITCASE[®] ramps allow durable, dependable solutions to a variety of mobility and accessibility needs. From the SUITCASE Singlefold to the Singlefold AS, Singlefold GF, and Trifold AS, you are sure to find a convenient, lightweight, portable, and easy-to-use ramp perfect for wheelchair and scooter use.

Visit SunriseMedical.ca for product specifications

Aluminum SUITCASE® Ramps

- Lifetime warranty
- Ideal for use at home or on the go
- Designed for wheelchairs and scooters
- Available in 3 styles, each with their own unique features
- Made in the USA



**SUITCASE®
Singlefold**

This portable ramp, with its single-fold design, sets up quickly and can be carried like a suitcase when not in use. Features a full-length hinge for additional strength and stability, and an extruded slip-resistant surface for superior traction. This durable ramp was designed for high traffic use and is simply the best in the industry. 800-pound weight capacity.



**SUITCASE®
Singlefold AS**

Can be used as one unit or easily separated into two lightweight sections, making it even more convenient to use, carry, and store. This ultra-strong ramp provides the safety and strength required while remaining light and portable. Features an applied, slip-resistant surface. 800-pound weight capacity.



**SUITCASE®
Trifold® AS**

Unique 3-fold design allows for ultimate portability when longer ramp lengths are needed. The ramp is designed to be used and carried as one unit, or can be quickly separated into two individual sections for easier transport and reduced handling weight. Features an applied, slip-resistant surface. 800-pound weight capacity.

Quick Specs

		Singlefold	Singlefold AS	Trifold AS
2'	Usable Size	2'x30"	2'x29-1/2"	—
	Folded Size	2'x15"x4"	2'x15"x4-1/2"	—
	Weight (lbs.)*	11	6 ea. (2 sections)	—
	Application	Curbs, 1 step	Curbs, 1 step	—
3'	Usable Size	3'x30"	3'x29-1/2"	—
	Folded Size	3'x15"x4"	3'x15"x4-1/2"	—
	Weight (lbs.)*	17	8 ea. (2 sections)	—
	Application	Curbs, 1 step	Curbs, 1 step	—
4'	Usable Size	4'x30"	4'x29-1/2"	—
	Folded Size	4'x15"x4"	4'x15"x4-1/2"	—
	Weight (lbs.)*	23	10 ea. (2 sections)	—
	Application	Curbs, 1-2 steps	Curbs, 1-2 steps	—
5'	Usable Size	5'x30"	5'x29-1/2"	5'x29-1/2"
	Folded Size	5'x15"x4"	5'x15"x4-1/2"	30-3/4"x15"x11-1/2"
	Weight (lbs.)*	29	13 ea. (2 sections)	16 ea. (2 sections)
	Application	Curbs, 1-2 steps, some minivans	Curbs, 1-2 steps, some minivans	Curbs, 1-2 steps, some minivans
6'	Usable Size	6'x30"	6'x29-1/2"	6'x29-1/2"
	Folded Size	6'x15"x4"	6'x15"x4-1/2"	36-3/4"x15"x11-1/2"
	Weight (lbs.)*	34	15 ea. (2 sections)	19 ea. (2 sections)
	Application	1-3 steps, some minivans & SUVs	1-3 steps, some minivans & SUVs	1-3 steps, some minivans & SUVs
7'	Usable Size	7'x30"	—	7'x29-1/2"
	Folded Size	7'x15"x4"	—	42-3/4"x15"x11-1/2"
	Weight (lbs.)*	41	—	22 ea. (2 sections)
	Application	1-3 steps, some minivans & SUVs	—	1-3 steps, some minivans & SUVs
8'	Usable Size	8'x30"	—	8'x29-1/2"
	Folded Size	8'x15"x4"	—	48-3/4"x15"x11-1/2"
	Weight (lbs.)*	47	—	25 ea. (2 sections)
	Application	2-3 steps, some minivans & SUVs	—	2-3 steps, some minivans & SUVs
10'	Usable Size	—	—	10'x29-1/2"
	Folded Size	—	—	60-3/4"x15"x11-1/2"
	Weight (lbs.)*	—	—	30 ea. (2 sections)
	Application	—	—	2-4 steps, some minivans & SUVs

Top Lip Extension (TLE)



Our optional TLE increases the reach of EZ-ACCESS SUITCASE ramps by extending the top lip from 3" to 9", allowing the ramp to clear the rear bumper of most SUVs and vans.

Be Prepared

Like a trusty umbrella or roadside emergency kit, your vehicle should always have a small wheelchair ramp tucked away for emergencies. A 2-ft. SUITCASE® Singlefold ramp is the ideal choice and its compact size allows for easy, take-anywhere portability. The lightweight ramp makes a great traveling companion because it fits perfectly in the trunk of a car, under the seat of a van, or is easily checked at the airport.

Call today for more information

mktcanada@sunmed.com | 1.800.263.3390 | www.SunriseMedical.ca

* Weights are approximate and may vary. Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.

Graphite Fiber SUITCASE® Ramps

- Lightest ramp on the market
- Lifetime warranty
- Ideal for use at home or on the go
- Designed for wheelchairs and scooters



SUITCASE® Singlefold GF

This ultra-lightweight ramp is exceptionally strong, thanks to its unique material composition consisting of glass-reinforced graphite fiber.

Available in 5 sizes, the SUITCASE Graphite Fiber ramp features a singlefold design and is the ideal choice for active users and caregivers.

Additional features include a full-length folding joint for additional strength, a slip-resistant coating for superior traction, and high visibility lines on the outer edges for added safety. 660-pound weight capacity.

For transportation on the back of a wheelchair, a carry bag for the 29" ramp is available.

29"	Usable Size	29" x 29" x 1"
	Folded Size	29" x 14.5" x 2"
	Weight (lbs.)*	7.75
	Application	Curbs, 1 step
35"	Usable Size	34.75" x 31" x 1"
	Folded Size	34.75" x 15.5" x 2"
	Weight (lbs.)*	9
	Application	Curbs, 1 step
51"	Usable Size	50.5" x 30" x 2"
	Folded Size	50.5" x 15.5" x 4.25"
	Weight (lbs.)*	13.25
	Application	Curbs, 1-2 steps
67"	Usable Size	66.5" x 30" x 2"
	Folded Size	66.5" x 15.5" x 4.25"
	Weight (lbs.)*	16.5
	Application	Curbs, 1-2 steps, some minivans
82"	Usable Size	81.75" x 30" x 2"
	Folded Size	81.75" x 15.5" x 4.25"
	Weight (lbs.)*	21
	Application	1-3 steps, some minivans & SUVs



mktcanada@sunmed.com | 1.800.263.3390 | www.SunriseMedical.ca

* Weights are approximate and may vary. Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.
Text and images © Homecare Products, Inc. All rights reserved.