




Gel "E" Skin Protection 3" Gel/Foam Wheelchair Cushion

Item #	Product Dimensions	Product HCPCS	UOM
14891	16" (W) X 16" (D) X 3" (H)	E2603	1/ea
14886	18" (W) X 16" (D) X 3" (H)	E2603	1/ea
14878	18" (W) X 18" (D) X 3" (H)	E2603	1/ea
14894	20" (W) X 16" (D) X 3" (H)	E2603	1/ea
14902	22" (W) X 18" (D) X 3" (H)	E2604	1/ea
14903	24" (W) X 18" (D) X 3" (H)	E2604	1/ea

- Pressure redistributing cushion designed to assist in the prevention, treatment and management of pressure injuries.
- Viscous gel bladder is encased in polyurethane foam shell and provides optimal pressure redistribution, support and comfort.
- Cover top is a urethane-coated nylon that is low-shear, vapor permeable and water-resistant.
- Cover base is waterproof vinyl for durability.
- Carry handle included for convenient transport.
- Meets flammability standard TB117.

specifications

Cover Material	Fluid Resistant Nylon Top Waterproof Vinyl Bottom
Weight Capacity	Up to 20" 275 lbs. 22" and up 500 lbs.
Warranty	18 Months Limited