



Groupe
Medicus

The Invacare Top End Pro Basketball Wheelchair can be quickly and easily adjusted in just 15 minutes with the quick release rear seat height adjustment system. Other adjustments include front seat height, back upholstery and footrest positioning, and center-of-gravity.

Adjustable platform footrests allows feet to be tucked where needed
Toe stop and strap keep feet in place
Available in tall (adult) or short (pediatric) frame
Welded 4-bend wing helps with the fast breaks
High-performance sports wheels with high pressure tires

Product Weight Capacity: 250 lb.