

GEMINO 20



GEMINO 20 • GEMINO 20 M • GEMINO 20 S

Rollator with unique functionality and flexibility

GEMINO 20 GEMINO 20 M GEMINO 20 S

Gemino 20 is available in three different sizes, so you can always be sure of finding a Gemino model that suits you and your needs.



IV pole



Cane holder



Light



Tray



Umbrella



Slow-down brake



Bag for oxygen bottle



Bag with lid and shoulder strap

Safety, support and stability

The unique and very stable design ensures that you can walk securely in all settings. When you need a rest, Gemino has a comfortable and convenient seat with back strap.

Comfortable, correct position for walking

It is easy to adjust Gemino to the correct height. The brake and pushing handles have been carefully designed to be comfortable and easy to use. A memory feature ensures handle heights can always return to the proper setting.

Safe and easy to use

High-quality integrated reflectors make Gemino clearly visible to road users. The unique brake block design ensures that the brakes stop the wheels very effectively. The design of the built-in curb climber makes it very easy to negotiate obstacles such as curbs and door sills. A unique side deflector makes it easier to navigate doorways, table legs and other furniture, without the risk of the back wheels getting caught.

Gemino – design

The items with which we surround ourselves help define how we are perceived. The elegant and sleek design of Gemino appeals to both young and old alike.

Specifications	Gemino 20	Gemino 20 M	Gemino 20 S
Total height	30.7 - 39.4 "	27.2 - 34.6"	25.6 – 30.3"
Total length	25.6 "	25.6"	25.6"
Total width	23.6 "	23.6"	21.7"

Distance between			
Pushing handles	18"	18"	16"
Seat height	24.4"	21.7"	18.9"
Total weight	16.1 lbs	15.9 lbs	15 lbs
Max. user weight	330 lbs	287 lbs	275 lbs
Max. weight in bag	11 lbs	11 lbs	11 lbs
Recommended user height	4'11" - 6'7"	4'5" - 5'7"	4'1" - 5'5"

Gemino is available in a number of different models:
Gemino 30 – a more advanced version of the rollator
Gemino 30 walker – with support for the forearms



The rollator folds up in one easy step and locks automatically into the folded position