

Platinum Collection Deluxe Padded Bath Seats

- Comfortable, cushioned seat and backrest are water-tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- Available with or without support arms
- Comes assembled with tool-free arm and back attachments
- Maximum Weight Capacity 300 lb EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0245
- · Lifetime limited warranty

7944A without support arms 7945A with support arms

1	ea
1	ea

<u>Specifications</u>	<u>7944A</u>	<u>7945A</u>	
Overall Width	20 ½"	23 ¾"	
Overall Depth	21"	21"	
Seat Width	15 ¾"	15 ¾"	
Seat Depth	16"	16"	
Minimum Seat Height	17 ½"	17 ½"	
Maximum Seat Height	22 ½"	22 ½"	
Height Of Back Support Above Seat	16"	16"	
Width Between Arms	N/A	17 ¾"	
Maximum Weight Capacity 300 lb EVENLY DISTRIBUTED			



