

Maple Leaf WHEELCHAIR



The only limitation is your imagination

HDCommode



STANDARD
HEAVY DUTY
COMMUNE



www.mapleleafwheelchair.ca

HDCommode

The HD COMMODE is a heavy duty commode that can be built to accomodate three levels of weight capacity and is always designed to fit the base of the frame over a standard toilet, even on custom orders.



♦ **COMMODE BUCKET**

COM1

WEIGHT CAPACITY - 250 LBS.

Includes:

- ♦ Seat to floor height - 22"
- ♦ Seat width - 18"
- ♦ Swing away footrests with composite footplates
- ♦ Black urethane arm pads & moulded seat
- ♦ Flip back armrests
- ♦ Back height - 20"
- ♦ Black vinyl padded back rest
- ♦ Black moulded sanitary cover
- ♦ (4) - 5" all brake casters
- ♦ 6" x 10" x 14" large sanitary bucket

COM2

WEIGHT CAPACITY - 500 LBS.

Includes:

- ♦ All the features listed in COM1 Package
- ♦ Seat width - 20"
- ♦ Overall width - 22"
- ♦ Heavy duty construction

COM3

WEIGHT CAPACITY - 750 LBS.

Includes:

- ♦ All the features listed in COM1 Package
- ♦ Seat width - 22"
- ♦ Overall width - 24"
- ♦ Extra heavy duty construction

ACCESSORIES & OPTIONS

- ♦ Custom seat width - 16", 17", 24"- 30"
- ♦ Custom seat height - 16"- 21" and 23"- 26"
- ♦ Custom seat depth - 16", 17", 19"- 22"
- ♦ Custom back heights - 21"- 24"
- ♦ Weight capacity up to 1000 lbs
- ♦ Custom soft seat
- ♦ Elevating legrests
- ♦ 24" rear wheels

*see full list of accessories on the HD COMMODE order form

MODEL SPECIFICATIONS

- ♦ 5 year warranty on frame
- ♦ 2 year warranty on additional components

Dealer Information

6540 Tomken Road
Mississauga ON L5T 2E9
T 905.564.2250 F 905.564.5424
1.800.362.0873

www.mapleleafwheelchair.ca

