

# Firefly Scooot

Scooot is a 4-in-1 mobility rider. It lets children with special needs play, explore and participate in ways that suit their physical and cognitive abilities.

- ✓ Scooot has four configurations; Scoot, Ride, Pooosh and Crawl.
- ✓ Suitable for children aged approx. 2-6 years old with a max. user weight of 22kg.
- ✓ Adjustable footrest to suit child's leg length as they grow.



### 4 in 1 Mobility Rider

Four different configurations- Scoot, Ride, Pooosh and Crawl- for different abilities and different ways to play.



### Independence

Scooot provides an early experience of Independence for lots of kids with mobility challenges. In doing so they will have opportunities for



### Advanced Backrest

For added support in Scooot, the Advanced Backrest Cushion can be added. The laterals move up and down to give support where it's needed, or they can be



### Adjustable Footplate

Adjustable footrest to suit the child's leg-length as they grow.

## Size Guide

Suitable for children aged approx. 2-6 years old with a max. user weight of 22kg.

Size	4-in-1 Configuration
Maximum User Weight	22kg/48.5lbs
Seat base width	230mm/9"
Seat base length	400mm/15.5"
Hip to end of footplate when seated	520mm/20.5" (not extended) 630mm/24.5" (fully extended)
Backrest width	260mm/10.5"
Backrest height	300mm/12"
Wheel arch width	230mm/9"