



STANDARD MODEL FEATURES

- Support users up to 66" in height and up to 130 lb in weight
- Low profile sides for easy transfer of users
- Durable, machine-washable padded cover
- Memory foam back/seat padding reduces sores and fatigue
- Soft-padded headwings
- Latex and silicone free
- Crash tested to FMVSS #213

PLUS MODELS FEATURE

1. Swing-Away Padded Trunk Supports
2. Swing-Away Padded Hip Supports
3. Swing-Away Abductor



UPHOLSTERY



Convertible Pink
(CP)



Go Kart Green
(GG)



Rally Blue
(RB)



Roadster Red
(RR)



Speedway Gray
(SG)

ACCESSORIES



Incontinence Covers



3-inch Extensor
Thrust Wedge



Buckle Guard



Retainer Clip Guard

ORDER FORM

FAX TO 905-787-8550 OR SUBMIT ONLINE

Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____

SIZE	Standard	Standard with Seat Extender	Plus - Medium	Plus - Large	Plus-M with Seat Depth Extender	Plus-L with Seat Depth Extender
Speedway Gray 	CSS-2400S-SG	CSS-2400E-SG	CSS-2400M-SG	CSS-2400L-SG	CSS-2401M-SG	CSS-2401L-SG
Roaster Red 			CSS-2400M-RR	CSS-2400L-RR	CSS-2401M-RR	CSS-2401L-RR
Convertible Pink 			CSS-2400M-CP	CSS-2400L-CP	CSS-2401M-CP	CSS-2401L-CP
Go Kart Green 			CSS-2400M-GG	CSS-2400L-GG	CSS-2401M-GG	CSS-2401L-GG
Rally Blue 			CSS-2400M-RB			CSS-2401L-RB
SELECT TO ITEMIZE		CSS-2400E-XX-IT	CSS-2400M-XX-IT	CSS-2400L-XX-IT	CSS-2401M-XX-IT	CSS-2401L-XX-IT
Seat Depth Extender		 CSS-K2425			 CSS-K2425	 CSS-K2425
Abductor			 CSS-K2415	 CSS-K2415	 CSS-K2419	 CSS-K2419
Hip Bracket			 CSS-K2473H (10" - 14")	 CSS-K2471H (13" - 16")	 CSS-K2473H (10" - 14")	 CSS-K2471H (13" - 16")
Hip Pad			 CSS-K2481M (4" x 8")	 CSS-K2481L (4" x 12")	 CSS-K2481M (4" x 8")	 CSS-K2481L (4" x 12")
Trunk Bracket			 CSS-K2474T (6" - 12")	 CSS-K2473T (10" - 14")	 CSS-K2474T (6" - 12")	 CSS-K2473T (10" - 14")
Trunk Pad			 CSS-K2482M (3" x 5")	 CSS-K2482L (4" x 6")	 CSS-K2482M (3" x 5")	 CSS-K2482L (4" x 6")
Weight Capacity	25 -130 lb	25 -130 lb	25 -130 lb	25 -130 lb	25 -130 lb	25 -130 lb
Max. Height Capacity	66"	66"	66"	66"	66"	66"
Harness Height	11", 13.5", 16", 18.5", 20.5", 23.5" (min. 10")					
Shoulder Width	18"	18"	18"	18"	18"	18"
Back Height	28"	28"	28"	28"	28"	28"
Seat Depth	12"	16"	12"	12"	16"	16"
Inside Seat Width at Front	17"	17"	17"	17"	19"	19"
Inside Seat Width at Rear	15"	15"	15"	15"	15"	15"
Overall Width	21.5"	21.5"	21.5"	21.5"	21.5"	21.5"
Seat-to-Back Angle	98°	98°	98°	98°	98°	98°
Product Weight	27 lb	29 lb	29 lb	29 lb	33 lb	33 lb
ACCESSORIES						
Incontinence Cover Kit	CSS-A2405S	CSS-A2405S	CSS-A2405M	CSS-A2405L	CSS-A2405M	CSS-A2405L
Extensor Thrust Wedge (raises 3" on the front)	CSS-A2455	CSS-A2455EX	CSS-A2455	CSS-A2455	CSS-A2455EX	CSS-A2455EX
Tether Strap Kit	CSS-A2030	one kit included in every car seat				
Bolt-in Hardware Kit	CSS-A2031	for vehicle manufactured before 2002				
Retainer Clip Guard	CSZ-A2085					
Buckle Guard	ZZ-A2075					