

## Luxury scooters

Shoprider's Crown series provides feature rich luxury in personal mobility. Experience the outdoors in unprecedented style, comfort and performance. The Crown Series are the most popular scooters in the Shoprider Canada portfolio.

### PathFinder 888<sup>SL/SE</sup>

- > **Special Edition** Shoprider 4-wheel scooter
- > Delta tiller with full light package
- > 1.2HP motor with 11.5 inch tires
- > Full suspension, anchor hooks
- > 50AH batteries with 4A charger
- > Colours: red, black or silver



### TrailBlazer 889<sup>SL</sup>

- > Shoprider 4-wheel scooter
- > Delta tiller with full light package
- > 1.2HP motor with 12.5 inch tires
- > Full suspension, anchor hooks
- > 50AH batteries with 4A charger
- > Colours: red, blue or silver



### TrailBlazer SE 889<sup>SL/SE</sup>

- > **Special Edition** Shoprider 4-wheel scooter
- > Delta tiller with 2-speed (hi-lo)
- > 2HP motor with 13.5 inch tires
- > Full suspension, full light package
- > 50AH batteries with 4A charger
- > Colours: red, black or silver



#### Crown series

Specs	Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chgr	Driving dist.	Turning radius	Safe climb°
888SL/SE	158KG/350LBS	48" x 23" x 43"	100KG/220LBS	1.2HP	12KPH/7.45MPH	2 x 50AH	4 x 11.5"	4AMP/HR	40KM/25M	1500MM/59"	9°/15%
889SL	158KG/350LBS	49" x 24.5" x 44"	104KG/229LBS	1.2HP	12KPH/7.45MPH	2 x 50AH	4 x 12.5"	4AMP/HR	40KM/25M	1500MM/59"	9°/15%
889SL/SE	158KG/350LBS	49" x 24" x 45"	100KG/220LBS	2.0HP	15KPH/9.45MPH	2 x 50AH	4 x 13.5"	4AMP/HR	40KM/25M	1500MM/59"	9°/15%