

Bellavita Dive Bath Lift

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

The lightest and lowest bath lift on the market now helps users dive deeper and get cleaner!

- New removable backrest with 50° recline
- Larger seat with hygiene cutout
- Improved dual column design provides optimal stability for up to 300 pounds
- Seat height adjusts from 2.6" to 18.9"
- Easy-to-use floating remote control
- Installs easy in almost any bathtubno plumbing or electrical work required

PRODUCTS	Description
477400252	Bellavita Dive Bath Lift w/ Wh Cvr, 1/EA
SPECIFICATIONS	Description
Largeur du produit	27.6"
CA Prop 65	DEHP chemical
Couleur	White
Matériel	Plastic
Longueur du produit	40.2"
Poids du produit	21.8 lbs
Dimensions de la boîte au détail	15.6" x 17.3" X 14.4"
Poids de la boîte au détail	5.8 lbs
Garantie	5yr Limited
Capacité de poids	300 lbs
Hauteur du produit	26.6