



4pt Sling bar Weight Scale Molift Mover 205

Medium 1530110

Large 1530111



Variants

	Item number	Compatible with	Height (mm)	Length (mm)	Safe working load (SWL) (kg)	Width (mm)
Weight scale SWL 300 kg	1100235	Molift Nomad	-	-	300	-
Weight Scale SWL 300 kg	1840000	Molift Air 200/205/300, Molift Mover 180	-	-	300	-
Weight Scale SWL 400 kg	1431021	Molift Air 350/500	-	-	400	-
2pt Sling bar 350 (S) incl. Weight Scale	1430120	Molift Mover 205	334	-	350	350
2pt Sling bar 450 (M) incl. Weight Scale	1430121	Molift Mover 205	334	-	350	450
2pt Sling bar 550 (L) incl. Weight Scale	1430122	Molift Mover 205	334	-	350	550
4pt Sling bar 450 (M) incl. Weight Scale	1530110	Molift Mover 205/300	439	450	350	300
4pt Sling bar 550 (L) incl. Weight Scale	1530111	Molift Mover 205/300	471	550	350	350