



Molift Battery & Charger

The Molift Power Pac is NiMH 14,4 V - 2,2 Ah and is suitable for the Molift Quickraiser 1, 2, 205 and 2+ as well as the Molift Smart 150 and Molift Mover 180.

It is the lightest, most powerful hoist battery on the market. Weighing as little as 1 Kg (2.2Lbs), this system will lift a 20+ stone person off the floor and into his chair, more than 10 times between charges.



Mover 205
Mover 300

Variants

	Item number	Battery	Compatible with
Molift Battery Charger	1340100	-	Charges NiMH battery packs with item no 0541000 (14,4V) and 1341000 (26,4V)
Molift Battery	0541000	NiMH 14.4 V - 2.6 Ah	Molift Quick Raiser 1 & 2, Molift Smart 150, Molift Mover 180
Molift Battery	1341000	NiMH 26.4 V - 2.6 Ah	Molift Mover 205/300, Molift Partner 255