



TRANSITIONS[®]

THRESHOLD RAMPS, MATS, and PLATES

Safety starts at your front door.



Visit SunriseMedical.ca for product specifications

Anytime. Anywhere. Safety starts at your front door.

EZ-ACCESS® knows that a gentle slope can make a world of difference. TRANSITIONS® threshold ramps, mats, and plates allow for a smooth ground-to-sill transition over doorway thresholds throughout your home. Reduce preventable trip hazards and make your home more comfortable and safe for mobility users.



The TRANSITIONS® Advantage

- Sturdy, seam-free surface makes TRANSITIONS ramps the safe solution for wheelchair, scooter, walker, and cane users
- Attractive designs are available in multiple styles and materials to suit your needs and taste
- Slip-resistant and designed for years of continued use
- Quick and easy setup and removal - installs in minutes so there's little-to-no disruption for you and your family
- Ideal for doorways, sliding glass doors, and raised landings
- Portable, convenient, and durable
- Manufacturer's warranty included
- Made in the USA

TRANSITIONS®
Threshold Ramps, Mats, and Plates




TRANSITIONS® Angled Entry Plate

- Formed aluminum construction for simple, one-step threshold navigation
- Three lengths are available (8, 10, and 12 inches) to accommodate thresholds up to 2 inches high
- Features an applied, slip-resistant surface for greater traction, even when wet
- 32 inches wide
- 700-pound weight capacity



TRANSITIONS® Angled Entry Ramp



- Swivel feet accommodate varying heights, angles, and surface types
- Three lengths are available (12, 24, and 36 inches) to accommodate thresholds ranging from 1-3/8 to 5-7/8 inches high
- Manufactured using high-strength aluminum that will not rust or corrode
- Features an extruded, slip-resistant surface for superior traction
- 36 inches wide
- 700-pound weight capacity

TRANSITIONS® Modular Entry Ramp

- Modular design with interchangeable extensions can accommodate thresholds from 3/4 to 6-1/4 inches high
- Installation is easy with our hidden hinge ... simply add and remove extensions to accommodate your needs
- Practical for sliding glass doors because it does not interfere with the door track
- Made of aluminum and features an extruded, slip-resistant surface for superior traction
- 34 inches wide
- 700-pound weight capacity



TRANSITIONS® Angled Entry Mat



- This "welcome" mat - made using recycled rubber - is eco-friendly and can be used safely on nearly any surface
- Stands up well in harsh outdoor conditions
- Features beveled edges for cross-traffic access
- Available in 1-1/2 and 2-1/2 inch heights and can easily be trimmed or notched if needed
- Usable width of 36 inches (40-43 inches overall)
- 850-pound weight capacity
- Available in three colors: Black, Storm Grey, and Hazelnut Brown

Call today for more information. We're happy to help.

mktcanada@sunmed.com | 1.800.263.3390 | www.SunriseMedical.ca

About EZ-ACCESS

EZ-ACCESS is the world's leading manufacturer of high-quality ramp systems. We are a family-owned and operated company with over 30 years of experience in designing and manufacturing access solutions. Our products are American made and are known for their safety, durability, and workmanship.

In addition to TRANSITIONS® Threshold ramps and mats, we also manufacture the versatile PATHWAY® Modular Access System, PASSPORT® Vertical Platform Lift, and a variety of portable ramps and other mobility accessories for both commercial and residential use.

EZ-ACCESS can make any type of home, regardless of age or style, accessible. Our products provide flexible short- and long-term solutions for a variety of needs.

Through an extensive network of local authorized dealers, EZ-ACCESS serves the entire United States with professional products and quick installation.

For additional information, please contact us today or visit ezaccess.com.



mktcanada@sunmed.com | 1.800.263.3390 | www.SunriseMedical.ca